

RESISTANCE THERMOMETER
Measuring insert: Fixed.

Type:
RT-RC1-K

Sheet No.
2-15 V2.1
5751-E010818V2.1




Application:

- For measuring temperatures in the food industry and the pharmaceutical, cosmetic and chemical/ technical industries

Properties:

- Sensor: Pt100, 2xPt100, Pt1000 in acc. with IEC 60751
- Measuring insert: Interchangeable
- Mechanical and thermal stress in accordance with DIN 43772
- Process attachment: Flange
- Outer protective sheath and flange: Stainless acid-proof steel
- Stands media temperatures of up to max. 250°C
- Quick reaction time
- IP67
- Approved by: GOST/TRCU



MECHANICAL SPECIFICATIONS

Protective sheath: -----
Wnr.1.4404 (AISI316L)

Nominal length L1: -----
20 mm
35 mm
Special

Process attachment: -----
Turning nipple G1/8 BSP
Turning nipple G3/8 BSP

Surface area of parts with media contact: -
Mechanically polished <0,8 µm
Electro-polished <0,5

ELECTRICAL SPECIFICATIONS:

----Sensor type:
1xPt100
2xPt100
1xPt1000

----Number of conductors:
2-wire
3-wire
4-wire

----Temperature range min/max:
-50/+250°C
Special

----Tolerance in acc. with IEC 60751:
Type A DIN (i.e.±(0,15+0,002xTactual)°C)
Type B 1/1 DIN (i.e.±(0,3+0,005xTactual)°C)
Type B 1/3 DIN(i.e.±(0,1+0,0017xTactual)°C)
Type B 1/6 DIN(i.e.±(0,05+0,00083xTactual)°C)
Type B 1/10 DIN(i.e.±(0,03+0,0005xTactual)°C)

----Cable type:
Cable type SS (Silicone-Silicone)
Cable type SBS (Silicone-Braid-Silicone)
Cable type TBT (Teflon-Braid-Teflon)
Special

----Cable length L3 [mm]
2 / 4 / 6 / 10 m
Special

Link for further information: [Pt100 Tolerance](#)

RESISTANCE THERMOMETER
Measuring insert: Fixed

Type:
RT-RC1-K

Sheet No.
2-15 V2.1
5751-E010818V2.1



Calibration:

Temperature calibration are used to verify and certify the sensor to have the correct accuracy. We can do either: "In house" or "Accredited" calibration. Accredited is certified by 3.e part. Normally we do a calibration in 3 points.

Enhanced performance services:


Cold applications (below -50°C) will influence the material and the measurement. CRYO treatment is needed to ensure a correct and working sensor down to -196°C.

A sensor will always drift over time, especially when there are high temperature fluctuations.

With "Ageing treatment" we stabilize the sensor to ensure a minimum drift over time. The benefits are long term stability, more correct measurement and easier planning of calibration periods.

Documentation:

Please order the correct documentation when ordering the sensor.



SIGNAL PROCESSING

Enclosure

Marine Box (112x82x42mm) -----

ABS Box (82x80x56mm) -----

NONE -----

Programmable mounted transmitter:

- [5333A Uninsulated for RTD](#)
- [5333D EEX Uninsulated for RTD](#)
- [5332A Uninsulated for RTD](#)
- [5332D EEX Uninsulated for RTD](#)
- [5331A Galvanic Isolated RTD / TC](#)
- [5331D EEX Galvanic Isolated RTD / TC](#)
- [5335A Hart 5 Protocol Standard](#)
- [5335D Hart 5 Protocol CSA, FM, ATEX, IECEx](#)
- [5337A Hart 5 & 7 Protocol](#)
- [5337D Hart 5 & 7 Protocol CSA, FM, ATEX, IECEx](#)



Transmitter Type:			
4 mA =	C°	20 mA =	C°

- [5350A Profibus standard](#)
- [5350B Profibus ATEX, FM and CSA](#)

Link to further information:

[Transmitter Overview](#)

[Programmable rail mounted transmitter](#)

CALIBRATION

---Calibration:

In house (Span -33°C - +700°C)

Accredited – in laboratory (-196°C - +1200°C)

1.	Point	°C
2.	Point	°C
3.	Point	°C

More point on request

Enhanced performance services

-----Cryo treatment.

For temperature sensor under -50°C

-----Ageing:

For long term stability.
Secure minimum drift of sensor accuracy

-----Documentation

- Certificate: 3.1 Material
- Certificate of origin
- Certificate of conformity
- Certificate of GOST

Other on request